

MICRO PUMP

MP-22 series

Continuous Micro Flow without Pulsation
to save limited / scarce samples

Innovative Pump Mechanism
with accurate control of valve switching



MP-22 series Specifications

Driving method	Liner-drive pump + auto switching valve
Selectable flow rate	0.1 – 999.9 μ L/min Note: Controllable range of pulsation : 10~500 μ L
Pump head material	Stainless Steel
Maximum pressure	20 MPa (2900psi)
Flow accuracy	< 0.3%
Communication	RS232C、RS485
Dimensions	150 (H) \times 160 (W) \times 280 (D) mm
Weight	Approx. 6.7 kg

FLOM Corporation

5-32-10 Shinmachi, Ome, Tokyo 198-0024 Japan
Phone: +81-428-30-7454
Website: www.flom.co.jp
Email: info@flom.co.jp

GL Sciences B.V.
www.glsciences.eu
info@glsciences.eu