

Pyrolysis Preparation Kit for OPTIC

GL Sciences Pyrolysis Preparation Kit (2406-4602) has a lot of convenient tools for your sample preparation. The Kit consists out of:

Reusable Material Punches (Various Sizes): These punches come in sizes ranging from 0.5 mm to 2.5 mm. They are ideal for creating consistent and precise sample pieces.

Cutting mats

- Large Cutting Mat (6 x 8 inches)
- Small Cutting Mat (2.5 x 3 inches)
- These mats provide a stable surface for cutting and handling samples.




Syringe for Aspirating Sample Pieces: Use this syringe to extract or transfer small sample pieces with precision.

Scalpel: A sharp scalpel is essential for precise cutting of samples.

DMI Sample Preparation Block:


- Designed for placing samples in the liners
- Includes 10 individual pieces



Tweezers: Fine-tipped tweezers are useful for handling delicate samples.

DMI Single Liner Weighing Block:

- Specifically designed for 5 mm liners
- Ensures accurate weighing of samples



All items are conveniently stored in a black suitcase.
Feel free to reach out if you need any further assistance: info@glsciences.eu